

Next Generation SedLine[®] Brain Function Monitor

More Complete Data, Now with an Enhanced Patient State Index (PSi)



Next Generation SedLine features:


- > An enhanced signal processing engine, which improves performance of the Patient State Index (PSi)
- > A Multitaper Density Spectral Array (DSA), which may enhance visibility of EEG features

Improved Patient State Index (PSi)

Next Generation SedLine features an enhanced signal processing engine which provides an enhanced Patient Sate Index (PSi), a processed EEG parameter related to the effect of anesthetic agents.

Expert Scoring¹

EEG experts scored the improvement in PSi performance between SedLine PSi and Next Generation SedLine PSi.




Experts found an overall 17% improvement in Next Generation PSi performance.

Next Generation SedLine Brain Function Monitoring helps clinicians monitor the state of the brain under anesthesia with bilateral data acquisition and processing of four leads of electroencephalogram (EEG) signals.




Parallel Signal Processing Engines

Next Generation SedLine utilizes Masimo's Parallel Signal Processing Engines to extract a clearer EEG signal for computing PSi.




The case below demonstrates Next Generation SedLine's improvement to PSi.¹




EMG is a common confounding factor that can interfere with EEG signals used in brain function monitoring.²

Adaptive Signal Processing with Band-Independent Features

When computing PSi, Next Generation SedLine uses adaptive signal processing with band-independent features to search for EEG features across many frequency bands.



The case below demonstrates Next Generation SedLine's improvement to PSi.¹




Power across all frequency bands decreases with age.³

Multitaper Density Spectral Array (DSA)


Next Generation SedLine offers clinicians the flexibility of choosing to display either an enhanced Multitaper Density Spectral Array (DSA) or a standard Hanning DSA. The DSA contains left and right spectrograms representing the power of the EEG on both sides of the brain.

When using a Multitaper DSA, EEG data are transformed into the frequency domain, which may provide a better display of EEG features.



Next Generation SedLine on Root®

The Next Generation SedLine module easily plugs into the Root patient monitoring platform via Masimo Open Connect® (MOC-9®) ports. Root's customizable, easily-interpretable display offers multiple views of brain monitoring information expanding visibility in the operating room and intensive care unit.



◀ **Next Generation SedLine** can be used simultaneously with O3® Regional Oximetry on the Root platform for a more complete picture of the brain.

RD SedLine™ EEG Sensor

- > Four active EEG leads collect data in the frontal lobe
- > Soft foam pads improve patient comfort
- > Allows simultaneous application of SedLine and O3 Regional Oximetry sensors



SedLine Specifications

PHYSICAL CHARACTERISTICS	ENVIRONMENTAL
Module Physical Dimensions Width 1.3 in (3.3 cm) Length 4.0 in (10.2 cm) Thickness 0.8 in (2.0 cm)	Module Operating Conditions Operating Temperature 41-104°F (5-40°C) Operational Humidity 15-95%, non-condensing Module Storage Conditions Storage Temperature -40-158°F (-40-70°C) Storage Humidity 15-95%, non-condensing Exposure to Pressure 500-1060 mbar

Sensor Specifications

Patient Weight > 30 kg	Ground Electrode CB
Application Site Forehead	Reference Electrode CT
Active Channels 4	Duration of Use Maximum of 24 hours
Active Electrodes L1, L2, R1, and R2	Latex Content Does not contain natural rubber latex

¹ Retrospective analysis of clinical data on file. ² Purdon P et al. *Brit J of Anaesth*. 10.1093 46-57. ³ Lobo, Francisco A., and Stefan Schraag. Limitations of anaesthesia depth monitoring. *Current Opinion in Anesthesiology*. 24, no. 6 (2011): 657-664.

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician. See instructions for use for full prescribing information, including indications, contraindications, warnings, and precautions.

Masimo U.S.
Tel: 1 877 4 Masimo
info-america@masimo.com

Masimo International
Tel: +41 32 720 1111
info-international@masimo.com

