



VACUUM REGULATORS



Turning the page on Vacuum Regulator products

VACUUM REGULATORS



VR-C3U2-F2A

ANALOG

- Analog gauge designed for clarity and appearance
- Operating ranges: 0-160 mmHg, 0-300 mmHg, and 0-760mmHg
- Upgrade to Intermittent or Digital in minutes
- Same housing for both Continuous and Intermittent regulator
- Removable backplate and gasket for ease of cleaning

DIGITAL

Includes all of the analog features plus...

- Patented state of the art Digital Display Technology
- Digital readout with no moving parts provides longer product life
- Accuracy within (+/-1%)
- Easy to read display



VR-C3U2-F2D

AVAILABLE MODELS

CONTINUOUS 2 MODE*



VR-C2U2-F2D

CONTINUOUS 3 MODE*



VR-C3U2-F2A

CONTINUOUS HIGH VAC*



VR-CHU2-F2D

CONTINUOUS/INTERMITTENT* (COMBO)



VR-CIU2-F2A

PEDIATRIC 2 MODE*



VR-P2U2-F2D

PEDIATRIC CONTINUOUS/INTERMITTENT* (COMBO)



VR-PIU2-F2A

**Options available:*

Analog Gauge



Selection of Colors




Digital Gauge




UPGRADES

ANALOG TO DIGITAL




CONTINUOUS TO INTERMITTENT



EASY TO SERVICE

- Can be sterilized internally with cold disinfectant
- Removable back plate for ease of cleaning
- 5 easy to replace modular components
- Service kits and parts are available



COLOR OPTIONS

CUSTOMIZE WITH COLOR

- Six additional colors to choose from, white is standard
- Different colors for ease of inventory management
- Colors are aesthetically pleasing



CUSTOMIZING EXAMPLE:



VR-C3U2-F2AB

ICU Blue

DIGITAL MODELS

CONTINUOUS MODELS

C2 - Continuous 2 Mode (0-300mmHg)
C3 - Continuous 3 Mode (0-300mmHg)
CH - Continuous High VAC (0-760mmHg)
P2 - Pediatric Continuous (0-160mmHg)



VR-CHU2-F2D

INLET-TO WALL

OUTLET-TO PATIENT

								
	No Fitting "F2"	Ohmeda "OM"	Chemetron "CM"	Puritan Bennett "PM"	DISS Handtite "DH"	DISS Nut "DN"	Oxequip "XM"	Medstar "MM"
	VR-C2U2-F2D	VR-C2U2-OMD	VR-C2U2-CMD	VR-C2U2-PMD	VR-C2U2-DHD	VR-C2U2-DND	VR-C2U2-XMD	VR-C2U2-MMD
	VR-C3U2-F2D	VR-C3U2-OMD	VR-C3U2-CMD	VR-C3U2-PMD	VR-C3U2-DHD	VR-C3U2-DND	VR-C3U2-XMD	VR-C3U2-MMD
	VR-CHU2-F2D	VR-CHU2-OMD	VR-CHU2-CMD	VR-CHU2-PMD	VR-CHU2-DHD	VR-CHU2-DND	VR-CHU2-XMD	VR-CHU2-MMD
	VR-P2U2-F2D	VR-P2U2-OMD	VR-P2U2-CMD	VR-P2U2-PMD	VR-P2U2-DHD	VR-P2U2-DND	VR-P2U2-XMD	VR-P2U2-MMD
	VR-C2UT-F2D	VR-C2UT-OMD	VR-C2UT-CMD	VR-C2UT-PMD	VR-C2UT-DHD	VR-C2UT-DND	VR-C2UT-XMD	VR-C2UT-MMD
	VR-C3UT-F2D	VR-C3UT-OMD	VR-C3UT-CMD	VR-C3UT-PMD	VR-C3UT-DHD	VR-C3UT-DND	VR-C3UT-XMD	VR-C3UT-MMD
	VR-CHUT-F2D	VR-CHUT-OMD	VR-CHUT-CMD	VR-CHUT-PMD	VR-CHUT-DHD	VR-CHUT-DND	VR-CHUT-XMD	VR-CHUT-MMD
	VR-P2UT-F2D	VR-P2UT-OMD	VR-P2UT-CMD	VR-P2UT-PMD	VR-P2UT-DHD	VR-P2UT-DND	VR-P2UT-XMD	VR-P2UT-MMD
	VR-C2UD-F2D	VR-C2UD-OMD	VR-C2UD-CMD	VR-C2UD-PMD	VR-C2UD-DHD	VR-C2UD-DND	VR-C2UD-XMD	VR-C2UD-MMD
	VR-C3UD-F2D	VR-C3UD-OMD	VR-C3UD-CMD	VR-C3UD-PMD	VR-C3UD-DHD	VR-C3UD-DND	VR-C3UD-XMD	VR-C3UD-MMD
	VR-CHUD-F2D	VR-CHUD-OMD	VR-CHUD-CMD	VR-CHUD-PMD	VR-CHUD-DHD	VR-CHUD-DND	VR-CHUD-XMD	VR-CHUD-MMD
	VR-P2UD-F2D	VR-P2UD-OMD	VR-P2UD-CMD	VR-P2UD-PMD	VR-P2UD-DHD	VR-P2UD-DND	VR-P2UD-XMD	VR-P2UD-MMD

CONTINUOUS/INTERMITTENT MODELS

CI - Continuous Intermittent - COMBO (0-300mmHg)
PI - Pediatric Continuous Intermittent - COMBO (0-160mmHg)



VR-PIU2-F2D

INLET-TO WALL

OUTLET-TO PATIENT

								
	No Fitting "F2"	Ohmeda "OM"	Chemetron "CM"	Puritan Bennett "PM"	DISS Handtite "DH"	DISS Nut "DN"	Oxequip "XM"	Medstar "MM"
	VR-CIU2-F2D	VR-CIU2-OMD	VR-CIU2-CMD	VR-CIU2-PMD	VR-CIU2-DHD	VR-CIU2-DND	VR-CIU2-XMD	VR-CIU2-MMD
	VR-PIU2-F2D	VR-PIU2-OMD	VR-PIU2-CMD	VR-PIU2-PMD	VR-PIU2-DHD	VR-PIU2-DND	VR-PIU2-XMD	VR-PIU2-MMD
	VR-CIUT-F2D	VR-CIUT-OMD	VR-CIUT-CMD	VR-CIUT-PMD	VR-CIUT-DHD	VR-CIUT-DND	VR-CIUT-XMD	VR-CIUT-MMD
	VR-PIUT-F2D	VR-PIUT-OMD	VR-PIUT-CMD	VR-PIUT-PMD	VR-PIUT-DHD	VR-PIUT-DND	VR-PIUT-XMD	VR-PIUT-MMD
	VR-CIUD-F2D	VR-CIUD-OMD	VR-CIUD-CMD	VR-CIUD-PMD	VR-CIUD-DHD	VR-CIUD-DND	VR-CIUD-XMD	VR-CIUD-MMD
	VR-PIUD-F2D	VR-PIUD-OMD	VR-PIUD-CMD	VR-PIUD-PMD	VR-PIUD-DHD	VR-PIUD-DND	VR-PIUD-XMD	VR-PIUD-MMD

COLOR SUFFIX

Your choice of colors by adding any of the following suffixes to the part number.
 White is standard if not otherwise specified.



Y - Yellow
 R - Red
 B - Blue
 G - Green
 P - Purple
 E - Grey

Example:
 VR-C3U2-F2DP



CONTINUOUS MODELS

C2 - Continuous 2 Mode (0-300mmHg)
C3 - Continuous 3 Mode (0-300mmHg)
CH - Continuous High VAC (0-760mmHg)
P2 - Pediatric Continuous (0-160mmHg)



VR-C3U2-F2A

INLET-TO WALL

OUTLET-TO PATIENT		INLET-TO WALL							
		No Fitting "F2"	Ohmeda "OM"	Chemetron "CM"	Puritan Bennett "PM"	DISS Handtite "DH"	DISS Nut "DN"	Oxequip "XM"	Medstar "MM"
		VR-C2U2-F2A VR-C3U2-F2A VR-CHU2-F2A VR-P2U2-F2A	VR-C2U2-OMA VR-C3U2-OMA VR-CHU2-OMA VR-P2U2-OMA	VR-C2U2-CMA VR-C3U2-CMA VR-CHU2-CMA VR-P2U2-CMA	VR-C2U2-PMA VR-C3U2-PMA VR-CHU2-PMA VR-P2U2-PMA	VR-C2U2-DHA VR-C3U2-DHA VR-CHU2-DHA VR-P2U2-DHA	VR-C2U2-DNA VR-C3U2-DNA VR-CHU2-DNA VR-P2U2-DNA	VR-C2U2-XMA VR-C3U2-XMA VR-CHU2-XMA VR-P2U2-XMA	VR-C2U2-MMA VR-C3U2-MMA VR-CHU2-MMA VR-P2U2-MMA
	Tubing Nipple "T"	VR-C2UT-F2A VR-C3UT-F2A VR-CHUT-F2A VR-P2UT-F2A	VR-C2UT-OMA VR-C3UT-OMA VR-CHUT-OMA VR-P2UT-OMA	VR-C2UT-CMA VR-C3UT-CMA VR-CHUT-CMA VR-P2UT-CMA	VR-C2UT-PMA VR-C3UT-PMA VR-CHUT-PMA VR-P2UT-PMA	VR-C2UT-DHA VR-C3UT-DHA VR-CHUT-DHA VR-P2UT-DHA	VR-C2UT-DNA VR-C3UT-DNA VR-CHUT-DNA VR-P2UT-DNA	VR-C2UT-XMA VR-C3UT-XMA VR-CHUT-XMA VR-P2UT-XMA	VR-C2UT-MMA VR-C3UT-MMA VR-CHUT-MMA VR-P2UT-MMA
	DISS "D"	VR-C2UD-F2A VR-C3UD-F2A VR-CHUD-F2A VR-P2UD-F2A	VR-C2UD-OMA VR-C3UD-OMA VR-CHUD-OMA VR-P2UD-OMA	VR-C2UD-CMA VR-C3UD-CMA VR-CHUD-CMA VR-P2UD-CMA	VR-C2UD-PMA VR-C3UD-PMA VR-CHUD-PMA VR-P2UD-PMA	VR-C2UD-DHA VR-C3UD-DHA VR-CHUD-DHA VR-P2UD-DHA	VR-C2UD-DNA VR-C3UD-DNA VR-CHUD-DNA VR-P2UD-DNA	VR-C2UD-XMA VR-C3UD-XMA VR-CHUD-XMA VR-P2UD-XMA	VR-C2UD-MMA VR-C3UD-MMA VR-CHUD-MMA VR-P2UD-MMA

CONTINUOUS/INTERMITTENT MODELS

CI - Continuous/Intermittent - COMBO (0-300mmHg)
PI - Pediatric Continuous/Intermittent - COMBO (0-160mmHg)



VR-CIU2-F2A

INLET-TO WALL

OUTLET-TO PATIENT		INLET-TO WALL							
		No Fitting "F2"	Ohmeda "OM"	Chemetron "CM"	Puritan Bennett "PM"	DISS Handtite "DH"	DISS Nut "DN"	Oxequip "XM"	Medstar "MM"
		VR-CIU2-F2A VR-PIU2-F2A	VR-CIU2-OMA VR-PIU2-OMA	VR-CIU2-CMA VR-PIU2-CMA	VR-CIU2-PMA VR-PIU2-PMA	VR-CIU2-DHA VR-PIU2-DHA	VR-CIU2-DNA VR-PIU2-DNA	VR-CIU2-XMA VR-PIU2-XMA	VR-CIU2-MMA VR-PIU2-MMA
	Tubing Nipple "T"	VR-CIUT-F2A VR-PIUT-F2A	VR-CIUT-OMA VR-PIUT-OMA	VR-CIUT-CMA VR-PIUT-CMA	VR-CIUT-PMA VR-PIUT-PMA	VR-CIUT-DHA VR-PIUT-DHA	VR-CIUT-DNA VR-PIUT-DNA	VR-CIUT-XMA VR-PIUT-XMA	VR-CIUT-MMA VR-PIUT-MMA
	DISS "D"	VR-CIUD-F2A VR-PIUD-F2A	VR-CIUD-OMA VR-PIUD-OMA	VR-CIUD-CMA VR-PIUD-CMA	VR-CIUD-PMA VR-PIUD-PMA	VR-CIUD-DHA VR-PIUD-DHA	VR-CIUD-DNA VR-PIUD-DNA	VR-CIUD-XMA VR-PIUD-XMA	VR-CIUD-MMA VR-PIUD-MMA

COLOR SUFFIX

Your choice of colors by adding any of the following suffixes to the part number.
 White is standard if not otherwise specified.



Y - Yellow
 R - Red
 B - Blue
 G - Green
 P - Purple
 E - Grey

Example:
 VR-C3U2-F2AB




ACCESSORIES

INLET / TO WALL



OUTLET / TO PATIENT



GAUGES




CANISTERS & BRACKETS



ROLL STANDS



DISTRIBUTED BY:



Toll Free: 1-866-GO-AMVEX (1-866-462-6839) • Tel: 905-764-7736 • Fax: 905-764-7743 • website: www.amvex.com
121 Granton Drive, Unit 21, Richmond Hill, Ontario, L4B 3N4, Canada